iMedPub Journals http://www.imedpub.com/

Journal of Clinical & Experimental Orthopaedics ISSN 2471-8416

Vol.4 No.1:50

DOI: 10.4172/2471-8416.100050

Safety and Efficacy of a Fibrin on Bleeding Bone

Birender Balain¹, Niall Craig², Kanna Gnanalingham³, Sanjeev Madan⁴, Hemant Sharma⁵, Guy Wynne-Jones⁶, John Hoseth⁷, Lian Øsytein⁷, Vidar Punsvik⁷, Sukumar Sura⁸, Rodolfo Padua⁹ and TimoFloydthy ^{9*}

¹Department of Orthopaedic Surgery, Robert Jones and Agnes Hunt Orthopaedic Hospital, Oswestry, England, UK

Received date: December 04, 2017; Accepted date: January 06, 2018; Published date: January 10, 2018

Citation: Balain B, Craig N, Gnanalingham K, Madan S, Sharma H, et al. (2018) Safety and Efficacy of a Novel Fibrin Dressing on Bleeding Cancellous Bone. J Clin Exp Orthop Vol 4, No 1: 50.

Copyright: ©2018 Balain B, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Abstract

Background: Cancellous bone bleeding (CBB) is a significant source of blood loss in spine and pelvis surgery. Current hemostatic methods are either ineffective or interfere with bone healing. A novel fibrin dressing (NFD) controls arterial, solid organ, and CBB in various species. The purpose of this clinical study was to determine the safety and performance of the NFD at controlling CBB.

Methods: Forty patients were enrolled at 9 centers in the UK, Norway and India. In each country Ethics Committee approvals were obtained prior to the study. This was a Level IV single-arm study of patients who met the inclusion eligibility criteria. The follow-up period was 6 weeks. The primary endpoint was control of bleeding at 3 min. Secondary endpoints were control of bleeding at 6 min, as well as clinical and laboratory criteria. Anatomical sites (AS) included posterior spinal fusion (PSF), iliac crest graft (ICG), and pelvic osteotomy (PO). Intraoperatively, the surgeon assessed the severity of CBB (pulsatile, flowing, oozing or none) and applied the dressing for 3 min. Bleeding was then assessed again and if it was not controlled (oozing or none) a second dressing was applied for an additional 3 min. For bleeding not controlled after the second dressing application the surgeon used an alternative method. Patients were assessed at 24 h and 6 weeks clinically (for hematoma, allergic reaction, etc.) and with lab analysis, including CBC, PT, aPTT and INR.

Results: All 40 patients completed the 6 week follow up. Twenty-three patients had PSF, 14 ICG and 3 PO. The dressing was used on a second AS in 11 of these patients for a total of 51 AS tested. Of the 51 treated AS, bleeding was controlled at 3 min in 41/51 (80.4%) and 6 min in 48/51 (94.1%). Three patients required alternative treatment. Bleeding in all 3 of the PO patients, the most vigorous bleeding, was controlled with the dressing. Adverse events (AE) considered to be "possibly device-related" occurred in 4 patients (13.3%).

Conclusions: This Phase I clinical trial demonstrated the safety and efficacy of this novel fibrin dressing, SurgiClot in cancellous bone bleeding. A randomized, controlled FDA trial currently is underway.

Introduction

Cancellous bone bleeding (CBB) is a significant source of blood loss in orthopedic and other surgery involving bone. Large blood loss in orthopedic surgery has been associated with transfusion risk, sequelae of shock, coagulopathy, electrolyte abnormalities, disseminated intravascular coagulation, infection, hematoma and neurologic injury [1-8]. Current hemostatic methods either are ineffective or interfere with bone healing [9-14]. Commonly used xenograft gelatin-based hemostatic agents swell, which can cause neurologic injury, anaphylaxis or intravascular thrombosis [15-17]. Other agents utilizing gelatin or cellulose have been useful in vascular surgery [18].

²Woodend Hospital & BMI Albyn Hospital, Aberdeen Royal Infirmary, Aberdeen, Scotland, UK

³Department of Neurosurgery, Salford Royal Foundation Trust, Manchester, England, UK

⁴Department of Orthopaedic Surgery, Doncaster and Bassetlaw Foundation Trust, Doncaster, England, UK

⁵Department of Orthopaedic Surgery, Spire Hull and East Riding Hospital, Hull, England, UK

⁶Department of Neurosurgery, Spire Washington Hospital, Newcastle, England, UK

⁷Department of Orthopaedic Surgery, Helse Møre og Romsdal HF, Kristiansund, Norway

⁸Department of Neurosurgery, Star Hospitals, Hyderabad, India

⁹Scientific Advisory Board, St. Teresa Medical, Inc., Eagan, MN, USA

^{*}Corresponding author: Timothy Floyd, MD, FACS, Chief Scientific Officer, Chairman, Scientific Advisory Board, St. Teresa Medical, Inc., Eagan, MN, USA, Tel: 001 208 841 2646; Fax: 001 208 367 7507; E-mail: t.floyd@stteresamedical.com

Vol.4 No.1:50

A novel fibrin dressing (trade name SurgiClot) controls bleeding in a variety of tissues in preclinical animal studies, including arterial and bone bleeding [19-21]. We conducted a single-arm clinical study to determine its safety and efficacy on CBB in humans as a Phase I prelude to a Phase III controlled trial. This is a report of our findings.

Materials and Methods

This was a Level IV single-arm study of patients who met the inclusion eligibility criteria. Enrollment was consecutive and was conducted in 9 centers in the UK, Norway and India. The follow-up period was 6 weeks. Forty patients were enrolled in this study.

When confronted intraoperatively with significant cancellous bone bleeding, the surgeon identified the site as appropriate for testing the NFD. Severity of CBB was graded according to a standardized system developed by Renkens. Bleeding was graded as either pulsatile, flowing, oozing or none. The surgeon then applied the dressing directly to the bleeding site and applied slight manual pressure using a standard gauze sponge. After 3 min the gauze sponge was removed and the bleeding site was assessed for the presence of residual bleeding which was assessed using the same criteria. If bleeding persisted, a second dressing was applied for another 3 min and then the site was reassessed. If significant bleeding persisted the surgeon then used other means to control the bleeding and the dressing was considered a failure.

The primary endpoint was control of bleeding at 3 min. Secondary endpoints included control of bleeding at 6 min, as

well as clinical and laboratory criteria. Treatment sites included posterior spinal fusion (PSF), iliac crest graft (ICG), and pelvic osteotomy (PO).

Patients were assessed clinically at 24 h and 6 weeks (for hematoma, allergic reaction, etc.) and with lab analysis, including complete blood count (CBC), pro-time (PT), activated partial thromboplastin time (aPTT) and international normalized ratio (INR). Surgeons completed a questionnaire regarding handling properties of the dressing.

Results

Thirty-eight of 40 patients completed the 6 week follow-up. Thirty-two PSF treatment sites were tested in 23 patients. Sixteen ICG treatment sites were tested in 14 patients. Three PO treatment sites were tested in 3 patients. The dressing was used at 2 anatomic sites in 11 of the 40 patients for a total of 51 treatment sites.

Of the 51 treatment sites, bleeding was controlled at 3 min in 41 (80.4%) and 6 min in 48 (94.1%). Nine of the 51 treatment sites required a second dressing application (17.6%) and a surgeon deviated from protocol on one patient who required a third application (2.0%) which controlled the bleeding. The second dressing application achieved hemostasis in 7 of the 9 sites (77.8%). Two sites required alternative treatment, but in total 49 of 51 (96%) bleeding sites were controlled with one or two dressing applications.

Table 1. Analysis of pre-op and post-op PT, aPT and INR.

	PT (sec)		aPTT (sec)		INR	
	Baseline	6-Week Follow-up	Baseline	6-Week Follow-up	Baseline	6-Week Follow-up
n	34	32	39	37	40	38
Mean ± SD	11.7 ± 0.9	11.7 ± 0.9	30.6 ± 4.2	32.9 ± 7.7	1.0 ± 0.1	1.0 ± 0.1
Median	11.5	11.5	31.5	31.7	1.0	1.0
Range	10-14.4	10-15	18.5-39	25-73	0.9-1.22	0.9-1.2

Interestingly, 100% of the secondary treatments sites (11/11) achieved hemostasis with a single dressing application at 3 min. Bleeding in all 3 of the PO patients, the most vigorous bleeding, was controlled with the dressing.

We saw no difference between pre-op and post-op protime (PT), activated partial thromboplastin time (aPTT) or international normalized ratio (INR) (Table 1).

Adverse events (AE) considered by the Clinical Review Committee (CEC) to be "possibly device-related" occurred in 2 patients (10%). These included elevated ESR and eosinophilia in 1 patient, monocytosis in 1 patient, neutrophilia and monocytosis in 1 patient, and wound infection in 1 patient.

A further secondary endpoint was handling characteristics determined by the surgeons (Table 2).

Table 2. Characteristics determined by the surgeons.

Characteristics	All Subjects % (n/N)
How easy is product to use?	·
Very easy	72.5% (29/40)
Easy	27.5% (11/40)
No opinion	0.0% (0/40)
Difficult	0.0% (0/40)
Very difficult	0.0% (0/40)

Very easy	65.0% (26/40)
Easy	35.0% (14/40)
No opinion	0.0% (0/40)
Difficult	0.0% (0/40)
Very difficult	0.0% (0/40)
How well did material conform to tissue surfac	ce?
Material conformed very well to tissue	60.0% (24/40)
Material conformed well to tissue	27.5% (11/40)
No opinion	5.0% (2/40)
Material conformed poorly to tissue	7.5% (3/40)
Material conformed very poorly to tissue	0.0% (0/40)
How well did the product dissolve?	-
Material dissolved very well	80.0% (32/40)
Material dissolved well	17.5% (7/40)
No opinion	2.5% (1/40)
Material dissolved poorly	0.0% (0/40)
Material dissolved very poorly	0.0% (0/40)
How effective was SURGICLOT at controlling topical haemostatic agents?	bleeding compared to o
Much more effective	47.5% (19/40)
More effective	32.5% (13/40)
About as effective	20.0% (8/40)
Less effective	0.0% (0/40)
Much less effective	0.0% (0/40)

Discussion

This feasibility study successfully demonstrated the efficacy and safety of this novel fibrin dressing.

We have shown it to be superior to a standard U.S. Army dressing in achieving hemostasis in arterial injuries in hypothermic, coagulopathic swine and it easily controlled bleeding cancellous bone in goats [19,20].

The surgeons participating in the present study were naive to the dressing yet were able to achieve 3 min hemostasis with one dressing application 80.4% of the time. They were able to achieve 3 min hemostasis 100% of the time upon treatment of subsequent different anatomic sites. Therefore, the learning curve is very short. Surgeons universally graded the dressing as either "as easy to use" or "easier to use" than currently available hemostatic products.

We believe it will have wide application in all surgical procedures involving bone and will help mitigate some of the complications and costs associated with blood loss in surgery where bone is cut, such as orthopedic surgery, spine surgery, cardiovascular surgery and head/neck surgery.

Since the dextran component quickly dissolves this dressing has the added advantage of leaving no residual foreign material. Pre-clinical safety studies have demonstrated inflammatory response comparable to controls with minimal foreign body reaction.

Conclusions

SurgiClot® is efficacious and safe at controlling CBB in spine and pelvis procedures (PSF, ICG, PO). A global randomized, controlled FDA trial of the dressing on bleeding cancellous bone currently is underway.

Ethics Committee Approvals

- 1. **Norway:** REK (Regionale Komiteer for Medicinisk of Helsefaglig Forskningsetikk) 2015-1382-16.
- 2. **United Kingdom:** WoSRES (West of Scotland Research Ethics Service) 14/WS/1153.
- 3. India: STAR Hospitals Institutional Ethics Committee SH/EC/AX2/SOP02/V1.

References

- Irita K (2011) Risk and crisis management in intraoperative hemorrhage: Human factors in hemorrhagic critical events. Korean J Anesthesiol 60: 151–160.
- Browne JA, Adib F, Brown TE, Novicoff WM (2013) Transfusion rates are increasing following total hip arthroplasty: risk factors and outcomes. J Arthroplasty 28: 34–37.
- Hu SS (2004) Blood loss in adult spinal surgery. Eur Spine J 13: S3– S5.
- Cha CW, Deible C, Muzzonigro T, Lopez-Plaza I, Vogt M, et al. (2002) Allogeneic transfusion requirements after autologous donations in posterior lumbar surgeries. Spine 27: 99–104.
- Möller H, Hedlund R (2000) Instrumented and noninstrumented posterolateral fusion in adult spondylolisthesis--a prospective randomized study: part 2. Spine 25: 1716–1721.
- Atwal NS, Bedi G, Lankester BJ, Campbell D, Gargan MF (2008) Management of blood loss in periacetabular osteotomy. Hip Int 18: 95–100.
- Pulido LF, Babis GC, Trousdale RT (2008) Rate and risk factors for blood transfusion in patients undergoing periacetabular osteotomy. J Surg Orthop Adv 17: 185–187.
- Dimitriou R, Mataliotakis GI, Angoules AG, Kanakaris NK, Giannoudis PV (2011) Complications following autologous bone graft harvesting from the iliac crest and using the RIA: a systematic review. Injury 42: S3–S15.
- Lavigne M, Boddu Siva Rama KR, Doyon J, Vendittoli PA (2008) Bone wax granuloma after femoral neck osteoplasty. Can J Surg 51: 58–60.
- Schonauer C, Tessitore E, Barbagallo G, Albanese V, Moraci A (2004) The use of local agents: bone wax, gelatin, collagen, oxidized cellulose. Eur Spine J 13: S89–S96.
- Nelson DR, Buxton TB, Luu QN, Rissing JP (1990) The promotional effect of bone wax on experimental Staphylococcus aureus osteomyelitis. J Thorac Cardiovasc Surg 99: 977–980.

Vol.4 No.1:50

- Nooh N, Abdullah WA, Grawish Mel-A, Ramalingam S, Javed F, et al. (2014) The effects of surgicel and bone wax hemostatic agents on bone healing: An experimental study. Indian J Orthop 48: 319– 325.
- Brodbelt AR, Miles JB, Foy PM, Broome JC (2002) Intraspinal oxidised cellulose (Surgicel) causing delayed paraplegia after thoracotomy--a report of three cases. Ann R Coll Surg Engl 84: 97– 99.
- 14. Ibarrola JL, Bjorenson JE, Austin BP, Gerstein H (1985) Osseous reactions to three hemostatic agents. J Endod 11: 75–83.
- Buchowski JM, Bridwell KH, Lenke LG, Good CR (2009) Epidural spinal cord compression with neurologic deficit associated with intrapedicular application of hemostatic gelatin matrix during pedicle screw insertion. Spine 34: 473–477.
- Spencer HT, Hsu JT, McDonald DR, Karlin LI (2012) Intraoperative anaphylaxis to gelatin in topical hemostatic agents during anterior spinal fusion: a case report. Spine J 12: e1–e6.
- Steinestel K, Geiger A, Naraghi R, Kunz U, Danz B, et al. (2013) Fatal thromboembolism to the left pulmonary artery by

- locally applied hemostatic matrix after surgical removal of spinal schwannoma: a case report. Hum Pathol 44: 294–298.
- Janczak D, Rucinski A, Rucinska Z, Janczak D, Dorobisz K, et al. (2013) Modern topical hemostatic agents—a breakthrough in vascular surgery. Polim Med 43: 141-146.
- Floyd CT, Rothwell SW, Martin R, Risdahl J, Olson CE (2012) A salmon thrombin-fibrinogen dressing controls hemorrhage in a swine model compared to standard koalin-coated gauze. J Spec Oper Med 12: 49–55.
- Floyd CT, Rothwell SW, Risdahl J, Martin R, Olson C, et al. (2012) Salmon thrombin-fibrinogen dressing allows greater survival and preserves distal blood flow compared to standard kaolin gauze in coagulopathic swine with a standardized lethal femoral artery injury. J Spec Oper Med 12: 16–26.
- Floyd CT, Padua RA, Olson CE (2017) Hemostasis and safety of a novel fibrin dressing versus standard gauze in bleeding cancellous bone in a caprine spine surgery model. Spine Deform 5: 310–313.